

MEDI TEKS



PRODUCT CATALOG

www.mediteks.com



CONTENTS

- Absorbent Gauze Bandage
- Conforming Light Elastic Bandage
- Med flex 3PLY Disposable Face Mask
- Hygiene Kit
- Cotton Wool
- Cotton Swabs
- Crepe Elastic Bandage (PBT)
- Compression Elastic Bandage (PBT)
- Coloured Compression Elastic Bandage
- Self Cohesive Bandage
- Cold&Hot Bandage
- Triangular Bandages
- Padding Bandages
- Plaster Of Paris Bandage (Gypsona)
- Absorbent Cotton Gauze
- Non-Sterile Gauze Compress
- Sterile Gauze Compress (Medical Pouch & Blister Packaging)
- Non-Sterile & Sterile Cotton Pad
- Non-Sterile & Sterile Abdominal (Laparotomy)Sponges
- First Aid Kit &First Aid Bag
- Instant Hot & Instant Cold Compress
- Vaseline Gauze Compress
- Re-Usable Hot & Cold Compress
- Re-Usable Hot & Cold Mask
- Re-Usable Hot & Cold Full Face Mask
- Re-Usable Hot & Cold Forehead Band
- Re-Usable Hot & Cold Chest Compress
- Re-Usable Hot & Cold Mickey
- Re-Usable Hot & Cold Knee Wrap
- Re-Usable Hot & Cold Wrist Wrap
- Re-Usable Hot & Cold Eye Pad
- Round Gauze Balls
- Sterile/Non Sterile Vaginal Tampons
- Sterile/Non Sterile Mikulicz Tampons



Established in 1999, MEDITEKS has made it its priority to excel continuously without compromising the company's basic values whilst complying with the rules of Quality Management. It has been the ultimate goal of the company to maintain increase total customer satisfaction.

The use of high-quality raw materials, coupled with advanced technology that is updated regularly, allows control at each and every stage of production MEDITEKS, whose product and service quality is approved both at home and abroad by official certificates, will continue to supply the best quality products at the most competitive prices to its valued clients around the world.



Minister of Health of the Republic of Turkey, Mr. Mehmet MÜEZZİNOĞLU and the Governor of Edirne, Mr. Hasan DURUER visiting our second, state of the art factory situated in the city of Edirne, on the European border and getting information from the senior management.

Minister of Health of the Republic of Turkey, Mr. Mehmet MÜEZZİNOĞLU and the Governor of Edirne, Mr. Hasan DURUER touring our Edirne factory and examining the MEDITEKS products with company managers.



The Mediflex Absorbent Gauze Bandage is produced from 100 %, first grade non-allergenic yarn. They are both much softer and more elastic so there are no loose threads. When used in open injuries in sterilized conditions, the Mediflex Absorbent Gauze Bandage minimizes the risk of infection. The Mediflex Absorbent Gauze Bandage is hydrophile (absorbent) so it is highly absorbent in case of bleeding.

PRODUCT CODE	PRODUCT SIZE	PACKAGING TYPE	NUMBER of PIECES / BOX	BOX DIMENSION
121202	5cm x 4m	Cellophane - Plastic Bag	3.000	76x46x52
121203	5cm x 5m	Cellophane - Plastic Bag	3.300	76x46x52
121204	5cm x 10m	Cellophane - Plastic Bag	2.000	76x46x52
121205	5cm x 100m	Cellophane - Plastic Bag	160	76x46x52
1212338	6cm x 5m	Cellophane - Plastic Bag	2.500	76x46x52
121207	7cm x 4m	Cellophane - Plastic Bag	2.400	76x46x52
121208	7cm x 5m	Cellophane - Plastic Bag	2.200	76x46x52
121209	7cm x 10m	Cellophane - Plastic Bag	1.000	76x46x52
121210	7cm x 100m	Cellophane - Plastic Bag	120	76x46x52
121577	8cm x 5m	Cellophane - Plastic Bag	2.000	76x46x52
121212	10cm x 4m	Cellophane - Plastic Bag	2.200	76x46x52
121213	10cm x 5m	Cellophane - Plastic Bag	700	50x50x30
121214	10cm x 10m	Cellophane - Plastic Bag	400	50x50x30
121215	10cm x 100m	Cellophane - Plastic Bag	120	76x46x52
1212339	12cm x 5m	Cellophane - Plastic Bag	1.600	76x46x52
121217	15cm x 4m	Cellophane - Plastic Bag	1.500	76x46x52
121218	15cm x 5m	Cellophane - Plastic Bag	600	50x50x30
121219	15cm x 10m	Cellophane - Plastic Bag	360	50x50x30
121220	15cm x 100m	Cellophane - Plastic Bag	40	50x50x30
121222	20cm x 4m	Cellophane - Plastic Bag	1.000	76x46x52
121223	20cm x 5m	Cellophane - Plastic Bag	1.100	76x46x52
121224	20cm x 10m	Cellophane - Plastic Bag	550	76x46x52
121225	20cm x 100m	Cellophane - Plastic Bag	60	76x46x52

*Please be informed that the product is not cutted from gauze. Side edges are not open it is specially weaved for secure use.



The Mediflex Conforming Light Elastic Bandage is produced by using viscose / cotton in weft thread and polyamide in warp thread. Both sides are woven by weaving machines so there are no loose threads. If required, this bandage can also be sterilized. Since the products woven by the weaving machines are both much softer and more elastic, the Mediflex Conforming Elastic Bandage can be wound in any desired length after being processed in the steam-boiler and the drying machine. The Mediflex Conforming Elastic Bandage has approximately 100 % elasticity. The Mediflex Conforming Elastic Bandage can be used in all areas where light-level support is required

PRODUCT CODE	PRODUCT SIZE	PACKAGING TYPE	NUMBER of PIECES / BOX	BOX DIMENSION
121227	6cm x 4m	Cellophane - Plastic Bag	2.200	76x46x52
121228	8cm x 4m	Cellophane - Plastic Bag	1.800	76x46x52
121229	10cm x 4m	Cellophane - Plastic Bag	1.700	76x46x52
121231	15cm x 4m	Cellophane - Plastic Bag	1.200	76x46x52
121232	20cm x 4m	Cellophane - Plastic Bag	800	76x46x52



- Fully Ultrasonic
- Standard Ear Loops
- 3 Fold Structure Design
- Face Mask With Filter Fabric
- Helps Prevent Diseases and Infections
- Filters Dust, Small Particular, Mist and Other Pollutants
- Soft and Flexible Nose Bridge, Can be Adjusted to the Face
- High Filtering Ability, Extremely Low Resistance, Easy to Breathe
- EN 14683:2019, EN 1041+A1:2014, EN 15223-1:2016, EN ISO 10933-1:2009/AC:2010, EN ISO 10993-5:2009 Test Approved

PRODUCT CODE	PRODUCT SIZE	PACKAGING TYPE	NUMBER of PIECES / BOX	BOX DIMENSION
1212437	3 ply Disposable Face Mask (meltblown)	50 pcs / box	2.500	51x41x51
1212438	3 ply Disposable Face Mask (spunbound)	50 pcs / box	2.500	51x41x51
1212135	3 ply Disposable Face Mask	1 pc / pack	-	-
1212136	3 ply Disposable Face Mask	5 pcs / pack	-	-
1212137	3 ply Disposable Face Mask	10 pcs / pack	-	-



Hygiene kit is ideal to keep in your handbag, car while travelling. Kit includes 5 pieces Meltblown face mask, 5 pieces cologne wipe, 50ml lemon cologne , 50ml anti-bacterial hand gel.

PRODUCT CODE	PRODUCT SIZE	PACKAGING TYPE	NUMBER of PIECES / BOX	BOX DIMENSION
1212526		12 pcs / pack		



Made from 100% Cotton fibers, bleached, absorbent, chemically neutral, odorless and without external impurities. It is designed for the needs of ambulances, hospitals and households.

PRODUCT CODE	PRODUCT SIZE	PACKAGING TYPE	NUMBER of PIECES / BOX	BOX DIMENSION
1212243	Cotton 1000 gr	PE Bag	24	760x210x1120
1212244	Cotton 500 gr	PE Bag	60	760x210x1120
1212245	Cotton 200 gr	PE Bag	75	600x990x300
1212246	Cotton 100 gr	PE Bag	150	600x990x300
1212247	Cotton 50 gr	PE Bag	300	600x990x300



Made from 100% Cotton fibers, bleached, absorbent, chemically neutral, odorless and without external impurities. It is designed for the needs of pharmacies, supermarkets, hospitals and households.

PRODUCT CODE	PRODUCT SIZE	PACKAGING TYPE	NUMBER of PIECES / BOX	BOX DIMENSION
1212248	40 pcs Round Make-Up Cleaning Pad	Pe Bag	50	410x510x400
1212249	100 pcs Round Make-Up Cleaning Pad	Pe Bag	50	410x510x400
1212250	100 pcs Make-Up Cleaning Pad	Pe Bag	48	270x340x350
1212251	120 pcs Make-Up Cleaning Pad	Pe Bag	48	340x470x300
1212252	150 pcs Make-Up Cleaning Pad	Pe Bag	48	340x470x300
1212253	100 pcs Cotton Swab	Plastic Box	144	285x415x513
1212254	200 pcs Cotton Swab	Plastic Box	144	505x505x350
1212255	300 pcs Cotton Swab	Plastic Box	144	395x595x505
1212256	200 pcs Square Cotton Swab	Plastic Box	144	465x550x300
1212257	60 pcs Square Box Cotton Swab for Babies	Plastic Box	144	465x550x300
1212258	90 pcs Square Box Cotton Swab for Babies	Plastic Box	144	460x535x430



The Mediflex Crepe Bandage is composed of polyester, cotton and latex. Since it is mostly first grade, 100 % cotton, it is non-allergenic and does not cause any itching. It can be used in all areas where medium-level support is required; such as muscular pains and bruises.

The Mediflex Crepe Bandage has 125 % elasticity and comes with two (2) metal or rubber clips – depending upon request - to fasten the bandage.

PRODUCT CODE	PRODUCT SIZE	PACKAGING TYPE	NUMBER of PIECES / BOX	BOX DIMENSION
121276	6cm x 150cm	Cellophane - Plastic Bag	1.700	76x46x52
121277	8cm x 150cm	Cellophane - Plastic Bag	1.500	76x46x52
121278	10cm x 150cm	Cellophane - Plastic Bag	1.200	76x46x52
121279	12cm x 150cm	Cellophane - Plastic Bag	900	76x46x52
121280	15cm x 150cm	Cellophane - Plastic Bag	850	76x46x52
121281	20cm x 150cm	Cellophane - Plastic Bag	550	76x46x52

* Elasticity: 125%

* Plastic and metal clips available upon request



Due to the high quality, natural fibers used, the Mediflex Compression Elastic Bandage is non-allergenic and eliminates the risk of itching. Its soft-strong texture permits blood circulation and thus it can easily be applied to all areas when required for long periods. The high quality raw materials used in the production of the Mediflex Compression Elastic Bandage make it washable and re-usable. The Mediflex Compression Elastic Bandage is ideal for use in varicosis, sportive traumas, soft-tissue tractions, dislocations and distortion bandaging, muscular and fiber captures. It has 150 % elasticity and comes with two (2) clips.

PRODUCT CODE	PRODUCT SIZE	PACKAGING TYPE	NUMBER of PIECES / BOX	BOX DIMENSION
121238	6cm x 1,50cm	Cellophane - Plastic Bag	1.100	76x46x52
121239	8cm x 1,50cm	Cellophane - Plastic Bag	900	76x46x52
121240	10cm x 1,50cm	Cellophane - Plastic Bag	800	76x46x52
121241	12cm x 1,50cm	Cellophane - Plastic Bag	550	76x46x52
121242	15cm x 1,50cm	Cellophane - Plastic Bag	750	76x46x52
121243	20cm x 1,50cm	Cellophane - Plastic Bag	400	76x46x52

* Elasticity: 150%

* Plastic and metal clips available upon request



Due to the high quality, natural fibers used, the Mediflex Coloured Compression Elastic Bandage is non-allergenic and eliminates the risk of itching. Its soft-strong texture permits blood circulation and thus it can easily be applied to all areas when required for long periods.

The Mediflex Coloured Compression Elastic Bandage can be used in the treatment of varicosis, orthopedic use, sportive traumas, soft-tissue tractions and surgeries. It has 150 % elasticity and comes with two (2) clips to fix the bandage.

PRODUCT CODE	PRODUCT SIZE	PACKAGING TYPE	NUMBER of PIECES / BOX	BOX DIMENSION
121250	6cm x 1,50cm	Cellophane - Plastic Bag	1.100	76x46x52
121282	8cm x 1,50cm	Cellophane - Plastic Bag	900	76x46x52
121283	10cm x 1,50cm	Cellophane - Plastic Bag	800	76x46x52
121284	12cm x 1,50cm	Cellophane - Plastic Bag	550	76x46x52
121285	15cm x 1,50cm	Cellophane - Plastic Bag	750	76x46x52
121286	20cm x 1,50cm	Cellophane - Plastic Bag	400	76x46x52



The Mediflex Self-Cohesive Bandage is self-cohesive so it is practical and easy-to-use. It is commonly used in accidents and as war first-aid. Its air-permeable and non-allergenic texture eliminates the risk of itching. Its soft-strong texture permits blood circulation and thus it can easily be applied to all areas when required for long periods. The Mediflex Self-Cohesive Bandage sticks only to itself. It is non-sticky to skin or hair. No pins or clips needed.

It is easy to dis-join. The Mediflex Self-Cohesive Bandage is mainly used in the treatment of varicosis, sportive traumas, strains, muscular captures and orthopedic use.

PRODUCT CODE	PRODUCT SIZE	PACKAGING TYPE	NUMBER of PIECES / BOX	BOX DIMENSION
1212178	2,5cm x 4m	Cellophane - Plastic Bag	500	29x30x29
121256	6cm x 4m	Cellophane - Plastic Bag	378	29x30x29
121267	8cm x 4m	Cellophane - Plastic Bag	252	29x36x29
121268	10cm x 4m	Cellophane - Plastic Bag	189	29x45x36
121290	15cm x 4m	Cellophane - Plastic Bag	126	30x33x60

* Only available in Brown.



The Mediflex Cold / Hot Bandage is a product that combines the advantages of Mediflex Cold / Hot Pack and Mediflex Elastic Bandage all in one. It is ideal for the treatment of sports injuries and all associated discomfort. Its special design molds to any part of the body perfectly so it is easily applicable. The Mediflex Cold / Hot Bandage is available in any desired length and colour. The Mediflex Cold / Hot Bandage mainly provides relief for sportive traumas, muscle pains and fiber captures. It is also highly recommended for usage in all areas where Cold / Hot Compress therapy is required.

PRODUCT CODE	PRODUCT SIZE	PACKAGING TYPE	NUMBER of PIECES / BOX	BOX DIMENSION
121251	15cm x 150cm	Boxed	25	62x25x28
121252	20cm x 150cm	Boxed	25	62x25x28

* Can be heated both in microwave oven and water



SPECIALY-DESIGNED PRODUCTS FOR PHARMACY DISTRIBUTION

- * Absorbent Gauze Bandage 10cmx1,5m (200 pcs / box)
- * Elastic Bandage (6cm, 8 cm, 10 cm, 12 cm ve 15 cm with box)
- * Sterile Gauze Swabs 30x80 (2 pcs 30x40) 150 pcs / box
- * Sterile Gauze Swabs 7,5x7,5 4 ply 5 pcs per bag (25 / 50 / 100 / box)

Re-Usable Cold / Hot Packs.(Please refer to pages 28-36)



The Mediflex Triangular Bandage is one of the most useful and versatile of medical bandages that can be utilized in an emergency. Amongst the many purposes of the Mediflex Triangular Bandage, the most common are arm-slings and tourniquets. They are also used to hold dressings in place splints together, to prevent infection and direct pressure.

PRODUCT CODE	PRODUCT SIZE	PACKAGING TYPE	NUMBER of PIECES / BOX	BOX DIMENSION
1212179	96cm x 135cm x 96cm	Cellophane - Plastic Bag	600	76x46x52



The Mediflex Padding Bandage is a light-weight bandage that offers protection of bony prominences and absorption of wound exudate. It is primarily used under other forms of compression bandaging and hard casts. The soft, texture of the Mediflex Padding Bandage bandage adheres to the cast and makes it possible to cut the cast easily to individual needs.

The highly extensible Mediflex Padding Bandage can be applied crease-free it does not slip and conforms to any body contour. The padding properties remain un-impaired over the whole period of wear. High water vapour and air permeability, as well as low absorption volume of the tightly-meshed bandage, guarantee a convincing wearer comfort. The high cotton content makes Mediflex Padding Bandage soft and very skin-friendly.

PRODUCT CODE	PRODUCT SIZE	PACKAGING TYPE	NUMBER of PIECES / BOX	BOX DIMENSION
1212175	10cm x 1,5m	Individual Cellophane	1.100	76x46x52
1212176	15cm x 1,5m	Individual Cellophane	800	76x46x52
1212177	20cm x 1,5m	Individual Cellophane	520	76x46x52
1212186	7,5cm x 2m	Individual Cellophane	800	76x46x52
1212180	10cm x 2m	Individual Cellophane	950	76x46x52
1212181	15cm x 2m	Individual Cellophane	550	76x46x52
1212182	20cm x 2m	Individual Cellophane	500	76x46x52
1212187	7,5cm x 2,7m	Individual Cellophane	700	76x46x52
1212183	10cm x 2,7m	Individual Cellophane	750	76x46x52
1212184	15cm x 2,7m	Individual Cellophane	510	76x46x52
1212185	20cm x 2,7m	Individual Cellophane	330	76x46x52
1212216	10cm x 3m	Individual Cellophane	700	76x46x52
1212217	15cm x 3m	Individual Cellophane	475	76x46x52
1212218	20cm x 3m	Individual Cellophane	360	76x46x52



The Mediflex plaster of paris bandage is a porous impregnated fast-setting plaster of paris bandage. On account of the special material structure made of a stable cotton carrier and the even application of the plaster, the Mediflex plaster of paris bandages can be moistened fast and offer excellent conformability and modelling properties.

The Mediflex plaster of paris bandages are available in following presentations: Bandages for casting with circular plaster of paris bandaging techniques; slabs as a splinting bandage providing additional strength applied on top of circular plaster of Paris bandages and large slabs for large plaster casts including plaster-of-Paris jackets as well as for extensive moulded casts. Setting time approx. 2-8 minutes.

PRODUCT CODE	PRODUCT SIZE	PACKAGING TYPE	NUMBER of PIECES / BOX	BOX DIMENSION
AS82125	5cm x 2cm	Individual Cellophane	150	38x27x25
AS82132	7,5cm x 2m	Individual Cellophane	125	38x27x25
1212138	10cm x 2m	Individual Cellophane	100	38x27x25
1212139	15cm x 2m	Individual Cellophane	100	38x27x25
1212140	20cm x 2m	Individual Cellophane	50	38x27x25
AS90486	5cm x 2,7cm	Individual Cellophane	125	38x27x25
AS82378	7,5cm x 2,7m	Individual Cellophane	110	38x27x25
1212161	10cm x 2,7m	Individual Cellophane	100	48x26x26
1212162	15cm x 2,7m	Individual Cellophane	50	37x26x26
1212163	20cm x 2,7m	Individual Cellophane	50	42x26x23
AS82552	5cm x 3cm	Individual Cellophane	125	42x26x23
AS82569	7,5cm x 3m	Individual Cellophane	100	42x26x23
1212171	10cm x 3m	Individual Cellophane	100	48x26x25
1212172	15cm x 3m	Individual Cellophane	50	37x26x26
1212173	20cm x 3m	Individual Cellophane	50	42x26x23



The Mediflex Absorbent Cotton Gauze is manufactured from 100% first grade non-allergenic absorbent cotton. It is white in colour and is highly absorbent which makes it widely used in hospitals as tampons and gauze swabs. The Mediflex Absorbent Cotton Gauze is the basic material for the production of Mediflex Gauze Swab, Mediflex Abdominal Sponge and Mediflex Sterile Cotton Pad.

PRODUCT CODE	PRODUCT SIZE	PACKAGING TYPE	NUMBER of PIECES / BOX	BOX DIMENSION
121287	90cm x 100m x 17 threads/cm ²	Bale	3.000	70x90x70
	90cm x 100m x 20 threads/cm ²	Bale	3.000	70x90x70
1212530	90cm x 350m	Roll		
1212531	90cm x 91m (36' x 100 yds) 4 ply	Roll		
1212532	90cm x 91m (36' x 100 yds) 4 ply x-ray	Roll		
1212533	90cm x 1000m	Roll		
1212534	90cm x 1.500m	Roll		

* One bale contains 30 rolls of 100 meters.



The Mediflex Non-Sterile Gauze Compress is produced from 100 %, first grade non-allergenic cotton yarn and thus is highly absorbent. It is extremely soft and gentle to the touch; yet strong and durable even when wet. The Mediflex Non-Sterile Gauze Compress is ideal for the general wound treatment, especially suitable in the primary treatment of dirty, infected and highly exudative wounds; as swabs and dressings in minor surgery in the outpatient department and on the ward. The cut edges are folded-in and the lack of any fraying minimizes the risk of any infection. The Mediflex Non-Sterile Gauze Compress is available in different plies and sizes.

PRODUCT CODE	PRODUCT SIZE	PACKAGING TYPE	NUMBER of PIECES / BOX	BOX DIMENSION
121269	5cm x 5cm 8 PLY	Cellophane - Plastic Bag	26.000	76x46x52
1212100	5cm x 5cm 12 PLY	Cellophane - Plastic Bag	20.000	76x46x52
121270	7,5cm x 7,5cm 8 PLY	Cellophane - Plastic Bag	9.000	65x46x50
1212164	7,5cm x 7,5cm 8 PLY X-RY	Cellophane - Plastic Bag	9.000	65x46x50
1212101	7,5cm x 7,5cm 12 PLY	Cellophane - Plastic Bag	10.000	76x46x52
1212165	7,5cm x 7,5cm 12 PLY X-RY	Cellophane - Plastic Bag	10.000	76x46x52
1212111	7,5cm x 7,5cm 16 PLY	Cellophane - Plastic Bag	9.000	76x46x52
1212112	7,5cm x 7,5cm 16 PLY X-RY	Cellophane - Plastic Bag	9.000	76x46x52
121271	10cm x 10cm 8 PLY	Cellophane - Plastic Bag	4.800	65x46x50
1212126	10cm x 10cm 8 PLY X-RY	Cellophane - Plastic Bag	4.800	65x46x50
1212102	10cm x 10cm 12 PLY	Cellophane - Plastic Bag	6.000	76x46x52
1212127	10cm x 10cm 12 PLY X-RY	Cellophane - Plastic Bag	6.000	76x46x52
1212113	10cm x 10cm 16 PLY	Cellophane - Plastic Bag	4.800	76x46x52
1212114	10cm x 10cm 16 PLY X-RY	Cellophane - Plastic Bag	4.800	76x46x52
1212166	8cm x 18cm 8 PLY	Cellophane - Plastic Bag	5.000	76x46x52
1212167	8cm x 18cm 8 PLY X-RY	Cellophane - Plastic Bag	5.000	76x46x52
1212115	8cm x 18cm 12 PLY	Cellophane - Plastic Bag	4.000	76x46x52
121272	8cm x 18cm 12 PLY X-RY	Cellophane - Plastic Bag	4.000	76x46x52

* All gauze swabs are folded with security folding method. There are no yarns coming out from sides.



The Mediflex Sterile Cotton Pad is manufactured from Mediflex Absorbent Cotton Gauze. The hydrophilic cotton, placed in the middle of the highly absorbent Mediflex Sterile Cotton Pad, makes it a multi-use cotton pad for medical dressings and cleansing procedures. The Mediflex Sterile Cotton Pad is widely used in obstetrics (as hygienic patient pads), in abdominal surgeries and as post-natal wound dressings.

PRODUCT CODE	PRODUCT SIZE	PACKAGING TYPE	NUMBER of PIECES / BOX	BOX DIMENSION
1212128	10cm x 10cm	Cellophane - Plastic Bag	1.920	76x46x52
121298	10cm x 15cm	Cellophane - Plastic Bag	1.360	76x46x52
1212123	10cm x 20cm	Cellophane - Plastic Bag	1.280	76x46x52
121299	12cm x 24cm	Cellophane - Plastic Bag	1.020	76x46x52

* On request it can be manufactured in a sterile package within different units



The Mediflex Sterile Gauze Compress is produced from 100 %, first grade non-allergenic cotton and thus is highly absorbent. It is conformable, tear-resistant and permeable. It is square-shaped and thus will adapt to wound surfaces of different sizes easily. The folded-in edges of the Mediflex Sterile Gauze Compress will help prevent fraying and minimize the risk of any infection. The Mediflex Sterile Gauze Compress is available in different plies and sizes.

PRODUCT CODE	PRODUCT SIZE	PACKAGING TYPE	NUMBER of PIECES / BOX	BOX DIMENSION
121258	30cm x 40cm	Medical Pouch	2.000	76x46x52
121259	30cm x 80cm	Medical Pouch	1.500	76x46x52
121260	5cm x 5cm (5 Pieces) 8 Ply	Medical Pouch	2.000	76x46x52
121261	7,5cm x 7,5cm (5 Pieces) 8 Ply	Medical Pouch	1.200	76x46x52
1212168	7,5cm x 7,5cm (25 Pieces) 4 Ply (5 Adet x 5 Package)	Boxed	324	76x46x52
1212169	7,5cm x 7,5cm (50 Pieces) 4 Ply (5 Pieces x 10 Pieces)	Boxed	180	76x46x52
1212170	7,5cm x 7,5cm (100 Pieces) 4 Ply (5 Pieces x 20 Pieces)	Boxed	90	76x46x52
121266	10cm x 10cm (5 Pieces) 8 Ply	Medical Pouch	1.000	76x46x52
1212124	8cm x 18cm (5 Pieces) 12 Ply	Medical Pouch	600	76x46x52
1212125	8cm x 18cm X-RY (5 Pieces) 12 Ply	Medical Pouch	600	76x46x52

* Different number of units can be packed and sterilized upon request.



PRODUCT CODE	PRODUCT SIZE	PACKAGING TYPE	NUMBER of PIECES / BOX	BOX DIMENSION
SGK83184	Sterile Gauze Compress 7,5cm x 7,5cm 10 pcs/pack	Blister Packaging	750	76x46x52
SGK83580	Sterile Gauze Compress 10cm x 10cm 8 Ply 10 pcs/pack	Blister Packaging	500	76x46x52
SGK83207	Sterile Gauze Compress 7,5cm x 7,5cm 8 Ply 10 pcs/pack	Blister Packaging	600	76x46x52
SBK94545	Sterile Lap Sponges 45x45 5pcs/pack	Blister Packaging	75	76x46x52

* Different number of units can be packed and sterilized upon request.



The Mediflex Abdominal Sponges are manufactured from 100 % surgical grade absorbent gauze. They are excellent for bleeding control and tissue protection during operations. The cut and folded-in edges of the Mediflex Abdominal (Laparotomy) Sponges will help prevent fraying and minimize the risk of any infection. The Mediflex Abdominal (Laparotomy) Sponges are available with or without a cloth anchoring loop and in different plies and sizes.

PRODUCT CODE	PRODUCT SIZE	PACKAGING TYPE	NUMBER of PIECES / BOX	BOX DIMENSION
1212116	20cm x 20cm 8 Ply X-ray	Cellophane-Plastic Bag	1.500	76x46x52
121273	30cm x 30cm 8 Ply X-ray	Cellophane-Plastic Bag	800	76x46x52
1212327	30cm x 30cm 12 Ply X-ray	Cellophane-Plastic Bag	600	76x46x52
121274	40cm x 40cm 8 Ply X-ray	Cellophane-Plastic Bag	400	76x46x52
NBK8440	40cm x 40cm 12 Ply X-ray	Cellophane-Plastic Bag	300	76x46x52
121275	45cm x 45cm 8 Ply X-ray	Cellophane-Plastic Bag	400	76x46x52
NBK71495	45cm x 45cm 12 Ply X-ray	Cellophane-Plastic Bag	300	76x46x52
1212129	20cm x 20cm 8 Ply	Cellophane-Plastic Bag	1.500	76x46x52
1212103	30cm x 30cm 8 Kat	Cellophane-Plastic Bag	800	76x46x52
1212326	30cm x 30cm 12 Ply	Cellophane-Plastic Bag	600	76x46x52
1212104	40cm x 40cm 8 Kat	Cellophane-Plastic Bag	400	76x46x52
NBK84433	40cm x 40cm 12 Ply	Cellophane-Plastic Bag	300	76x46x52
1212105	45cm x 45cm 8 Ply	Cellophane-Plastic Bag	400	76x46x52
NBK84419	45cm x 45cm 12 Kat	Cellophane-Plastic Bag	300	76x46x52

* Please be informed upon request different quantity in pouch can be sterilized.

* The product is produced with secure folding, side stitches prevent any yarn coming throughout the sides.



CONTENT

- 1 Mediflex Cold / Hot Pack 14x14 cm
- 1 pair of Scissors
- 1 pair of Examination Gloves
- First air Plasters (10 pieces)
- Absorbent Cotton gauzes Bandage 10X150 cm
- Elastic bandage 6 cm 1 pieces
- Conforming Elastic Bandage 10 cm x 2 m (1 pieces)
- Conforming Elastic Bandage 6 cm x 4 m (1 pieces)
- Vaseflex Sterile Paraffinned Gauze Swab (1 piece)
- Sterile Gauze Swab 30 cm x 40 cm (1 piece)

Mediflex First-Aid Kit is a well-stocked first-aid kit that is roomy, durable, easy to carry, and simple to open. Its plastic tackle box is lightweight, have handles, and offers a lot of space. Kept within easy reach, the well-endowed Mediflex First-Aid Kit will help you handle an emergency at a moment's notice.

PRODUCT CODE	PRODUCT SIZE	PACKAGING TYPE	NUMBER of PIECES / BOX	BOX DIMENSION
121244	-	-	88	76x46x52

First - Aid Bag



CONTENT

- Absorbent Cotton Gauze Bandage 5X150 cm (3 Pieces)
- Absorbent Cotton Gauze Bandage 7X150 cm (3 Pieces)
- Absorbent Cotton Gauze Bandage 10X150 cm (3 Pieces)
- Absorbent Cotton Gauze Bandage 10x4 m
- Airway
- Antiseptic Solvent
- Light Elastic Bandage 6 cm
- Elastic Bandage 6cm
- Sterile Gauze Swab 30 cm X 80 cm (30 x40 - 2 pieces)
- Sterile Gauze Swab 7.5 cm X 7.5 cm (10 pieces)
- Tourniquet
- Triangular Bandage
- Vaseflex (Sterile Paraffinned Gauze Swab 2 pieces)
- Instant Cold Pack
- Cotton Wool 50 gr
- 1 pair of Scissors
- First Aid Plaster (10 Pieces)
- Pair of Examination Gloves
- Scissors
- Safety Pin
- Notebook
- Pencil

PRODUCT CODE	PRODUCT SIZE	PACKAGING TYPE	NUMBER of PIECES / BOX	BOX DIMENSION
121245	-	-	35	76x46x52

Instant Cold Compress

MEDI—TEKS



Mediflex Instant Cold Pack is easily activated with a simple squeeze for emergency use where cold therapy is required and temporarily relieves minor pain and swelling for sprains, aches and sore joints. Mediflex Instant Cold Pack reaches 0 degrees Celsius (32 degrees Fahrenheit) within seconds after you squeeze the inner bag and release the water into the crystals in the outer bag. The cold lasts from about 20-25 minutes, depending upon the temperature of the area it is placed and the ambient temperature. Mediflex Instant Cold Pack is great for use within first aid kits and to keep in gym bags and stays at a lower temperature longer than other cold packs. Single use only.

PRODUCT CODE	PRODUCT SIZE	PACKAGING TYPE	NUMBER of PIECES / BOX	BOX DIMENSION
121253	13cm x 23cm	-	Boxed 25 Unboxed 44	29x27x17 63x25x30

MEDI—TEKS

Instant Hot Compress



Mediflex Instant Hot Pack is easily activated for emergency use where hot therapy is required and is commonly used to help in reducing post-operative bruising and relieving pain, soreness, tension, stress, and aches in the leg, chest, back and abdominal areas that can occur due to surgery or injury.

The Mediflex Instant Hot Pack is ready for use in less than a minute with a simple squeeze and twist of the bag. This odor-free, convenient and safe hot pack continues to deliver therapeutic heat for up to 20 minutes and stays at a warmer temperature longer than other hot packs. Single use only.

PRODUCT CODE	PRODUCT SIZE	PACKAGING TYPE	NUMBER of PIECES / BOX	BOX DIMENSION
121254	13cm x 23cm	-	Boxed 25 Unboxed 44	29x27x17 63x25x30

MEDI—TEKS

www.mediteks.com

26



VASETULLE is for use as a topical treatment only. It is used to cover the wounds, prevent it from contamination, keep it wet and accelerate healing process. The range of wounds on which VASETULLE may be used include:

- Minor burns and scalds.
- Skin loss wounds, lacerations and abrasions.
- Donor and recipient graft sites.

PRODUCT CODE	PRODUCT SIZE	PACKAGING TYPE	NUMBER of PIECES / BOX	BOX DIMENSION
P-021	Vasetulle 10cm x 10cm	Medical Paper	1.600 - 2.000	35x54.5x27.5cm
P-022	Vasetulle 30cm x 10cm	Medical Paper	400	35x36x28cm
P-011	Hexatulle 10cm x 10cm	Medical Paper	1.600 - 2.000	35x54.5x27.5cm
P-012	Hexatulle 30cm x 10cm	Medical Paper	400	35x36x28cm
P-031	Betatulle 10cm x 10cm	Medical Paper	1.600 - 2.000	35x54.5x27.5cm
P-032	Betatulle 30cm x 10cm	Medical Paper	400	35x36x28cm
P-051	Fucitulle 10cm x 10cm	Medical Paper	1.600 - 2.000	35x54.5x27.5cm



Mediflex Re-Usable Cold / Hot Pack is a dual-purpose compress that provides icy cool or soothing heat. It contains a state of the art gel which provides soothing warmth or cold for extended periods of time.

Mediflex Re-Usable Cold / Hot Pack is a versatile and practical compress that can be used on any part of the body. This is an economical product that can be used over and over again. In addition, it is available in many convenient sizes and a great number of fun and attractive colours.

Mediflex Re-Usable Cold / Hot Pack is used to reduce the pain in physical injuries, post medication traumas, sprains, second degree burns, headache and migraine, bruises, nose-bleed and toothache. It is also widely used to provide relief for tired and injured muscles, rheumatism, muscle- hardening, cricks, menstrual pain and slow-blood circulation. Ideal as a drug-free therapy for minor pain. Be sure to follow the instructions for use found on the box and the pack.

PRODUCT CODE	PRODUCT SIZE	PACKAGING TYPE	NUMBER of PIECES / BOX	BOX DIMENSION
121246	14cm x 14cm	Boxed	116	62x25x28
121247	12cm x 24cm	Boxed	60	62x25x28
121248	14cm x 28cm	Boxed	60	62x25x28
121257	26cm x 35cm	Boxed	10	62x25x28

*Can only be heated both in microwave oven and water.



Mediflex Re-Usable Cold / Hot Mask is a re-usable face and cosmetic mask with a velvety side that does no harm to skin. It is used on the face for both cold and hot therapy to relieve facial pain and to help in reducing post-operative swelling and bruising.

The compress may be used for a wide range of cosmetic procedures, including chemical peel, collagen injections, eyelid surgery, face lift, laser skin resurfacing, nose surgery, and other procedures. In addition, cold-compress is to relieve pain from headache and migraine, tired and swollen eyes, fever and hangover. Warm-compress is suggested when suffering from catarrh, influenza and sinus inflammation.

The gel contained provides a slow cooling or warming of the face in a natural, drug-free way.

PRODUCT CODE	PRODUCT SIZE	PACKAGING TYPE	NUMBER of PIECES / BOX	BOX DIMENSION
121288	-	Boxed	60	63x25x28

*Can only be heated both in microwave oven and water



Mediflex Re-Usable Cold / Hot Full-Face Mask is a dual-purpose, re-usable face and cosmetic mask that provides icy cool or soothing heat. It contains a state of the art gel which provides soothing warmth or cold for extended periods of time.

Mediflex Re-Usable Cold / Hot Full-Face Mask is widely used to help in reducing post-operative swelling and bruising from cosmetic procedures. It is ideal after face lift and botox treatment recommended by plastic surgeons following laser-surgery and microdermabrasion, for sunburns, migraines and headaches.

Mediflex Re-Usable Cold / Hot Full-Face Mask also stimulates circulation for a healthy look, clears clogged pores after steam application and soothed redness after facial peels. It can also be used after or during massage treatments.

PRODUCT CODE	PRODUCT SIZE	PACKAGING TYPE	NUMBER of PIECES / BOX	BOX DIMENSION
121291	-	Boxed	25	63x25x28

*Can only be heated both in microwave oven and water.



Mediflex Re-Usable Cold / Hot Forehead Band is a versatile hot or cold compress which contains a state of the art gel that provides soothing warmth or cold for extended periods of time. Its convenient and leakproof design conforms to eyes and distributes its weight evenly over the applied surface.

Mediflex Re-Usable Cold / Hot Forehead Band is widely used to provide relief from tired and puffy eyes, dark circles, swelling from injury or surgery, tension headache, nasal congestion, eye strain and hangover.

PRODUCT CODE	PRODUCT SIZE	PACKAGING TYPE	NUMBER of PIECES / BOX	BOX DIMENSION
121292	-	Boxed	90	62x25x28

*Can only be heated both in microwave oven and water.



The Mediflex Re-Usable Cold / Hot Breast Compress comes as a pair and may be used on the breast for both cold and hot therapy to help in reducing post-operative swelling and bruising. It also diminishes the pains that can occur while breast-feeding and after nursing.

The flexible non-cracking PVC thermo plastic membrane conforms gently to the breasts and the shoulder straps allow you to move while using the compress.

The Mediflex Re-Usable Cold / Hot Breast Compress can be used for a wide range of cosmetic uses, including breast augmentation, breast reduction, breast lift, breast reconstruction, medical tattooing, male breast reduction (gynecomastia), mastectomy, nipple reconstruction and in relief from pains associated with breast feeding and nursing.

PRODUCT CODE	PRODUCT SIZE	PACKAGING TYPE	NUMBER of PIECES / BOX	BOX DIMENSION
1212110	-	Boxed	40	62x25x28

*Can only be heated both in microwave oven and water



Mediflex Re-Usable Cold / Hot Mickey is ideal for emergency first aid needs at home or away and provides an excellent solution for those who prefer drug-free therapy for kids.

It is a dual-purpose, odor-free, re-usable compress covered with a material cover to protect children's delicate skin. Mediflex Re-Usable Cold / Hot Mickey is a versatile and practical compress that can be used on any part of the body. This is an economical product that can be used over and over again. In addition, it is available in a great number of fun and attractive colours.

Mediflex Re-Usable Cold / Hot Mickey is used to treat minor bumps, bruises, muscle aches, strains and sprains, tension headaches, toothaches, minor burns, scrapes, cuts, insect bites, respiratory tract discomfort, fever and infections.

PRODUCT CODE	PRODUCT SIZE	PACKAGING TYPE	NUMBER of PIECES / BOX	BOX DIMENSION
121289	-	Boxed	25	62x25x28

*Can only be heated both in microwave oven and water.



Mediflex Re-Usable Cold / Hot Knee Wrap is a dual-purpose, ergonomically-designed, re-usable and economic thermo-therapy wrap that provides icy cool or soothing heat. It contains an odor-free, state of the art gel which provides soothing warmth or cold for extended periods of time.

Its flexible non-cracking PVC thermo plastic membrane conforms gently to the knees and the adjustable straps allow you to move freely while using the wrap.

Mediflex Re-Usable Cold / Hot Knee Wrap is used to provide relief from sprains, pains, strains, bruises and common knee conditions including knee arthritis, ligament injuries, runner's knee and other associated knee problems.

PRODUCT CODE	PRODUCT SIZE	PACKAGING TYPE	NUMBER of PIECES / BOX	BOX DIMENSION
121293	18cm x 37cm	Boxed	25	62x25x28

*Can only be heated both in microwave oven and water.



Mediflex Re-Usable Cold / Hot Wrist Wrap is a dual-purpose, ergonomically-designed, re-usable thermo-therapy wrap that provides icy cool or soothing heat. It contains a state of the art gel which provides soothing warmth or cold for extended periods of time. Mediflex Re-Usable Cold / Hot Wrist Wrap is ideal for use during occupational, physical and sports therapy as well as simple day-to-day needs. It provides relief from post surgery discomfort, muscle strains and sprains effectively by offering both hot and cold treatment within one product. Mediflex Re-Usable Cold / Hot Wrist Wrap also helps reduce carpal tunnel syndrome, pain and muscle spasms.

The flexible non-cracking PVC thermo plastic membrane conforms gently to the wrist and on the pack adhesive straps allow you to move freely while using the wrap. Ideal as a drug-free therapy for minor pain. Be sure to follow the instructions for use found on the box and on the pack.

PRODUCT CODE	PRODUCT SIZE	PACKAGING TYPE	NUMBER of PIECES / BOX	BOX DIMENSION
121294	11cm x 24cm	Boxed	60	62x25x28

*Can only be heated both in microwave oven and water.



The Mediflex Re-Usable Cold / Hot Eye Mask can be used for cold or hot relief. Cold relieves dry, puffy eyes, allergies and sinus headaches. It helps take the sting out of a brow wax and helps reduce swelling after a cosmetic procedure. Heat helps to reduce nasal congestion and relieve soreness once the swelling is gone.

PRODUCT CODE	PRODUCT SIZE	PACKAGING TYPE	NUMBER of PIECES / BOX	BOX DIMENSION
1212131	-	Plastic Bag	100	29x27x17

*Can only be heated both in microwave oven and water.



The Mediflex Round Gauze Balls are produced from 100 %, first grade non-allergenic cotton yarn. They are hydrophilic and thus are mainly used to absorb blood and clean wounds during operations. There is no fraying so the risk of infection is minimized. The Mediflex Round Gauze Balls are also effective in the impediment of bacterial invasions.

PRODUCT CODE	PRODUCT SIZE	PACKAGING TYPE	NUMBER of PIECES / BOX	BOX DIMENSION
1212130	8cm x 8cm X-RY	Cellophane-Plastic Bag	100	-
1212116	10cm x 10cm X-RY	Cellophane-Plastic Bag	100	-
1212117	12cm x 12cm X-RY	Cellophane-Plastic Bag	100	-
1212118	15cm x 15cm X-RY	Cellophane-Plastic Bag	100	-
1212119	20cm x 20cm X-RY	Cellophane-Plastic Bag	100	-
1212120	25cm x 25cm X-RY	Cellophane-Plastic Bag	100	-
1212121	30cm x 30cm X-RY	Cellophane-Plastic Bag	100	-
1212122	40cm x 40cm X-RY	Cellophane-Plastic Bag	100	-



The Mediflex Vaginal Tampons are used to absorb blood from the body after vaginal operations. It is produced as sterile and non-sterile.

PRODUCT CODE	PRODUCT SIZE	PACKAGING TYPE	NUMBER of PIECES / BOX	BOX DIMENSION
1212383	-	Cellophane-Plastic Bag	400	76x46x52



Mikulicz Tampon is used mainly in brain surgeries and cystectomy. It is produced as sterile and non-sterile.

PRODUCT CODE	PRODUCT SIZE	PACKAGING TYPE	NUMBER of PIECES / BOX	BOX DIMENSION
1212384	-	Cellophane-Plastic Bag	1.760	76x46x52

Dotted lines for notes.

MEDİ — TEKS

MEDİTEKS SAĞLIK HİZMETLERİ TIBBİ MALZ. TEKS. ÜRÜN. SAN. ve DIŞ TİC. A.Ş.

Istanbul Headquarters

📍 Maslak Mah. Ahi Evran Cd. No:21, Polaris Plaza,
Kat 6, Office No: 33, 34398 Sarıyer / İstanbul - Türkiye

☎ +90 212 674 29 68

📠 +90 212 674 05 20

Edirne Factory

📍 Edirne Organize Sanayi Bölgesi 4. Cadde
No:5 Süloğlu / Edirne - Türkiye

☎ +90 284 316 23 81 - 82

📠 +90 284 316 23 84



mediteks@mediteks.com 🌐 www.mediteks.com